




1 Attach the connecting hoses to the PAD



2 Apply the PAD in the correct position



3 Attach the hoses to the unit




4 Fill the tank with liquid

## START SESSION WITH PREINSTALLED TREATMENTS



5 Select (1) quick set program  
Select (2) AIR compression settings



6 You can select 5 different pulsations of AIR  
compression settings pressing on button (2)




7 Select (3) preinstalled treatments  
(4) edit the treatment




8 If you click the button (4) you can set  
parameters temperature and duration

# START SESSION WITH QUICK START - NEW TREATMENT



If you click the button (1) or (3) you go to play screen set (5) change the temperature, set (6) the treatment duration and set (7) the pause duration **9**



(8) you can repeat treatment and (9) play the treatment **10**



Change temperature (10) even in execution **11**



You can select 3 different levels of AIR compression (11) **12**

HIGH (65-75 mmHg)

MEDIUM (45-55 mmHg)

LOW (15-20 mmHg)

NOTE: Tolerance +/- 10 mmHg in base of atmospheric pressure



Apply first the wrap to the body part art before using compression



Attention! Avoid bending the tubes at the wrap connectors



Advice: Between treatment sessions - keep the unit and pad connected.

## WARNINGS

1. Use only ZAMAR certified NON-TOX liquid 3LIQ010
2. Fill the liquid tank to the arrow marked on the label next to the liquid tank
3. In case of transport drain the liquid tank
4. Always be sure that the top plug is properly closed
5. Check always the seals on the plug and connectors
6. Do not tilt the unit upside down

1. Turn off the electrical switch before working on the unit
2. Ensure the plug load supports the unit power requirements (Amps)
3. Ensure the generator power supports the unit power requirement (kVA)
4. Do not share the outlet/circuits with other electrical power devices
5. Using an extension power cord can cause excessive voltage drops